



www.IetLtd.com Proudly serving laboratories worldwide since 1979

CALL +847.913.0777 for Refurbished & Certified Lab Equipment

Bruker Avance 300 AV300



BRUKER Magnetics

Low Loss Cryostats

Superconducting Magnets

phone: ++41 1 825 91 11

fax: ++41 1 825 92 15

e-mail: magnetics@bruker.ch

service@magnetics.bruker.ch

sales@magnetics.bruker.ch

SUPERCONDUCTING

ULTRASHIELD™

NMR MAGNET SYSTEM


Coil Number BZH 89/300/70A

Dewar Number D 232/54-3046



www.IetLtd.com Proudly serving laboratories worldwide since 1979

CALL +847.913.0777 for Refurbished & Certified Lab Equipment



UltraShield™ NMR Magnet System UltraShield™ 300 MHz/54 mm long hold

3 Superconducting NMR Magnet **300/70A**

3.1 Characteristic Data **89 - 300/70A**

Proton Frequency	300 MHz
Central Field	7.05 Tesla
Coil Inductance	14.4 Henry
Magnetic Energy	52 kJoule
Magnetic center from top flange *)	459 mm
Main Coil Heater Current	120 mA
Shim Coil Heater Current	150 mA

	Magnet-Test	System-Test	Customer Site
Magnet Current	A 83.7		83.8 8407
X-Shim Current	A +1.13		+2.0 +1.56
Y-Shim Current	A +7.66		+8.2 7.66
Z-Shim Current	A -7.27		-7.8 -8.0
XZ-Shim Current	A +0.96		+1.0 -0.72
YZ-Shim Current	A -0.14		-0.1 -2.04
XY-Shim Current	A -1.18		-1.1 -1.74
X ² -Y ² -Shim Current	A -1.81		-1.81
Z ² -Shim Current	A -5.93		-5.8 -5.1
Z ³ -Shim Current	A -3.32		-3.3 -3.38
Frequency change due to Z ² -Shim	kHz +5.49		
Magnetic center from top flange	mm 464		464
RT Shim System Angle **)	Deg 90°		90°
Visa	PMA		

Important: During charging Z and Z² shim heaters must be permanently ON

Remarks: Magnetsystem is NMR tested!

*) Approximate values (mechanical drawings) after cool down.
 **) Measured from the right-hand He stack to the cable input of the RT shim system.



UltraShield™ NMR Magnet System

UltraShield™ 300 MHz/54 mm long hold

3.3 Resistance Measurements

88-300/70A

Measurements at room temperature with the current lead mounted in the cryostat:					
from	A	Connector A			
to	L	Connector B	39,5	OHM	Main Heater
from	C	Connector A			
to	L	Connector B	4,0	OHM	Z Heater
from	E	Connector A			
to	L	Connector B	4,0	OHM	X Heater
from	F	Connector A			
to	L	Connector B	3,9	OHM	Y Heater
from	H	Connector A			
to	L	Connector B	3,9	OHM	XZ Heater
from	J	Connector A			
to	L	Connector B	3,8	OHM	YZ Heater
from	K	Connector A			
to	L	Connector B	3,8	OHM	XY Heater
from	L	Connector A			
to	L	Connector B	3,6	OHM	X ² -Y ² Heater
from	D	Connector A			
to	L	Connector B	3,8	OHM	Z ² Heater
from	K	Connector B			
to	L	Connector B	3,4	OHM	Z ³ Heater
from	A,B	Connector B			
to	D,E	Connector B	4,3	OHM	Shim Coils +/-
from	+	High Curr. Conn.			
to	H	Connector B	1,0	OHM	High Current to Sense +
from	+	High Curr. Conn.			
to	-	High Curr. Conn.	7,4	OHM	Main Coil
from	-	High Curr. Conn.			
to	J	Connector B	0,9	OHM	High Current to Sense -
from	H	Connector B			
to	J	Connector B	7,9	OHM	Sense + Sense -
from	A,B	Connector B			
to	L	Connector B	>30M	OHM	Shim Coil to Heater common
from	D,E	Connector B			
to	H	Connector B	>30M	OHM	Shim Coil to Maincoil
from	H	Connector B			
to	L	Connector B	>30M	OHM	Sense to Heater common
from		the Connectors			
to		the Ground	>30M	OHM	Insulation Magnet to Dewar
At room temperature with connection lead for the ACD* mounted in the cryostat:					
from	K	Connector ACD*			
to	J	Connector ACD*	111	OHM	Upper Temperature Sensor PT100
from	A	Connector ACD*			
to	B	Connector ACD*	110	OHM	Lower Temperature Sensor PT100