



IET | International Equipment Trading Ltd.
www.IetLtd.com | Proudly serving laboratories worldwide since 1979
 CALL +847.913.0777 for Refurbished & Certified Lab Equipment

Dual Stage Rotary Vane Pumps

Performance Specifications

Technical Specifications	Hz	Units	DS 102	DS 202	DS 302	DS 402	DS 602	DS 1002
Free air displacement	60	l/min (cfm)	114 (4)	192 (6.8)	285 (10)	410 (14.5)	605 (21.4)	1010 (35.7)
	50	l/min (m3/h)	95 (5.7)	160 (9.6)	237 (14.2)	342 (20.5)	504 (30.2)	842 (50.5)
Pumping speed	60	cfm	3.5	5.8	8.2	12.3	18.3	31
	50	m3/h	5	8.3	11.6	17.4	26	44
Ultimate partial pressure		mbar	10, -4	10, -4	10, -4	10, -4	10, -4	10, -4
Ultimate total pressure		mbar	2 x 10 ⁻³	2 x 10 ⁻³	2 x 10 ⁻³	2 x 10 ⁻³	2 x 10 ⁻³	2 x 10 ⁻³
Ultimate total pressure w/gas ballast		mbar	2 x 10 ⁻²	2 x 10 ⁻²	2 x 10 ⁻²	2 x 10 ⁻²	2 x 10 ⁻²	2 x 10 ⁻²
Water vapor tolerance		mbar	15	15	20	30	30	20
Water vapor capacity		g/h	60	100	160	350	550	650
Oil capacity min/max		l	0.2/0.5	0.25/0.6	0.25/0.6	0.5/1	0.5/1	2.3/3.1
Motor rating 1ph	50/60	kW	0.45/0.55	0.45/0.55	0.45/0.55	0.75/0.90	0.75/0.90	
Motor rating 3ph	50/60	kW				0.75/0.90	0.75/0.90	1.10/1.30
Nominal rotation speed	60	rpm	1800	1800	1800	1800	1800	1800
	50	rpm	1500	1500	1500	1500	1500	1500
Weight		kg	22	25	25	35	35	62
		lb	48	55	55	77	77	136
Inlet flange		DN	25KF	25KF	25KF	25KF	25KF	40KF
Exhaust flange		DN	25KF	25KF	25KF	25KF	25KF	40KF

Exhaust Filters - Oil Mist Eliminator

The exhaust and the oil mist eliminator retain the oil vapors which would otherwise be expelled into the atmosphere during pumpdown and gas ballast operation. Some models feature an oil return line to allow

IET - Used Lab Equipment - Refurbished Analytical Laboratory Instruments

condensed oil or fluid to return to the pump reservoir. Varian's new traps are designed to prevent the backstreaming of mechanical pump fluids. Copper and stainless steel gauze inserts are designed to reduce oil backstreaming. Molecular sieve inserts are available for applications where it is desirable to increase water vapor pumping speed while eliminating backstreaming above the trap.

Model	A	B	C	D	H	L	M	N	P	R	W
DS 102	30 (1.18)	30 (1.18)	35 (1.38)	3 (0.12)	212 (8.3)	430 (16.9)	245 (9.6)	104 (4.1)	105 (4.13)	165 (6.5)	132 (5.2)
DS 202	30 (1.18)	30 (1.18)	35 (1.38)	3 (0.12)	212 (8.3)	467 (18.4)	245 (9.6)	104 (4.1)	105 (4.13)	205 (8.1)	132 (5.2)
DS 302	30 (1.18)	30 (1.18)	35 (1.38)	3 (0.12)	212 (8.3)	467 (18.4)	245 (9.6)	104 (4.1)	105 (4.13)	205 (8.1)	132 (5.2)
DS 402	40 (1.57)	40 (1.57)	40 (1.57)	6 (0.24)	242 (9.5)	541 (21.3)	261 (10.3)	130 (5.1)	100 (3.9)	230 (9)	164 (6.5)
DS 602	40 (1.57)	40 (1.57)	40 (1.57)	6 (0.24)	242 (9.5)	541 (21.3)	261 (10.3)	130 (5.1)	100 (3.9)	230 (9)	164 (6.5)

Model	A	B	C	H	L	M	N	P	R	W	Z
DS 1002	65 (2.6)	55 (2.2)	95 (3.7)	335 (13.2)	640 (25.2)	355 (14)	185 (7.3)	115 (4.5)	300 (11.8)	230 (9)	350 (13.8)

Replacement Parts Info:

Dual Stage Rotary Vane Pumps	Part Number
DS 102 with 1 phase worldwide motor*	9499315
DS 102 with 1 phase worldwide motor* (UK version)	9499350
DS 202 with 1 phase worldwide motor*	9499320
DS 202 with 1 phase worldwide motor* (UK version)	9499351
DS 302 with 1 phase worldwide motor*	9499325
DS 302 with 1 phase worldwide motor* (UK version)	9499352
DS 402 with 1 phase worldwide motor*	9499330
DS 402 with 3 phase worldwide motor**	9499331
DS 602 with 1 phase worldwide motor*	9499335
DS 602 with 3 phase worldwide motor**	9499336
DS 1002 with 3 phase worldwide motor**	9499340

1 phase motors (100-120 / 200-240) V \pm 10%, 50/60 Hz. ** 3 phase motors (200-220 / 380-415) V \pm 10% at 50 Hz or (200-230 / 460) V \pm 10% at 60 Hz. All motors comply with CE and UL/CSA standards.

Minor Maintenance Kit	Part Number
(Contains all the valves, O-rings and seals to refurbish the pump to vacuum integrity)	
DS 102, DS 202, DS 302	9499370
DS 402, DS 602	9499371
DS 1002	9499372

Major Maintenance Kit	Part Number
(Includes all the items of the minor kit plus the vanes)	
DS 102	9499380
DS 202, DS 302	9499381
DS 402, DS 602	9499382
DS 1002	9499383

Oil and Accessories	Part Number
Rotary vane fluid, GP type, 1 liter	9499390
Rotary vane fluid, GP type, 1 liter (USA)*	K7516301
Rotary vane fluid, GP type, 2 liters	9499391
Rotary vane fluid, GP type, 1 gallon (USA)	K7516302
Oil mist eliminator (DS 102 - DS 602)	9499395
Oil mist replacement cartridge (pkg. of 2)	9499394
NW 25 oil exhaust filter	9499392
NW 40 oil exhaust filter	9499393
NW 25-NW 40 oil exhaust replacement cartridge	9499342
Oil return kit	9499376
European plug power cable 2 meters (DS 102 - DS 602) 1 ph	9499396
USA plug power cable 2 meters (DS 102 - DS 602) 1 ph	9499397
UK plug power cable 2 meters (DS 102 - DS 602) 1 ph	9499398
NW 25 centering ring Viton (DS 102 - DS 602)	KC25AV
NW 25 aluminum clamp (DS 102 - DS 602)	KQ25AWP
NW 40 centering ring Viton (DS 1002)	KC40AV

NW 40 aluminum clamp (DS 1002)

KQ40AWP

* When these pumps are used in Leak Detectors applications, we recommend the use of Rotary Vane Fluid (Elite-Z mechanical), P/N 695409005 as it features a lower vapor pressure.



WWW.IETLTD.COM



IET | International Equipment Trading Ltd.

www.IetLtd.com | *Proudly serving laboratories worldwide since 1979*

CALL +847.913.0777 for Refurbished & Certified Lab Equipment

IET - Used Lab Equipment - Refurbished Analytical Laboratory Instruments